# **fitvibe**smart<sup>®</sup>

*Your personal coach*

User Manual

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Mode d'emploi

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<b>English</b>	.....	<b>3</b>
<b>Deutsch</b>	.....	<b>19</b>
<b>French</b>	.....	<b>37</b>
<b>Português</b>	.....	<b>53</b>
<b>Español</b>	.....	<b>67</b>

# **fitvibe smart**<sup>®</sup>

*Your personal coach*

## User Manual

## **User Manual Fitvibe Smart**

### **Device for controlled Whole Body Vibration (cWBV)**

Manufacturer: Uniphy Elektromedizin GmbH & Co. KG  
Main office: Neuendorfstrasse 19b  
D-16761 HENNIGSDORF  
Telephone: +49 (0)33025044-0  
Fax: +49 (0)33025044-99  
E-mail: [info@fitvibe.com](mailto:info@fitvibe.com)  
Website: [www.fitvibesmart.com](http://www.fitvibesmart.com)  
[www.fitvibe.com](http://www.fitvibe.com)

### **© 2007, Uniphy Elektromedizin GmbH & Co. KG**

All rights reserved. Nothing from this publication or the fitvibe smart device may be copied, stored in an automated data file, or made public, in any form or in any way, be it electronically, mechanically, by photocopying, recordings or in any other way, without prior written permission from Uniphy Elektromedizin GmbH & Co. KG

Art. Code 331.837

Version 1.2  
June 2007

# Table of Contents

<b>1</b>	<b>Safety</b>	<b>7</b>
1.1	Purpose	7
1.2	Legal aspects of training	7
1.3	Safety instructions	7
1.3.1	General	7
1.3.2	Contra-Indications	7
1.3.3	Unexpected symptoms	8
1.3.4	Electrical Safety	8
1.3.5	Prevention of explosion	8
1.3.6	Liability	8
1.3.7	Conformity	8
1.3.8	Deactivation and disposal	9
<b>2</b>	<b>Components and Installation</b>	<b>9</b>
2.1	Receipt	9
2.2	Description of the control unit and accessory items	10
2.3	Assemblage	11
2.4	Transport & Storage	11
<b>3</b>	<b>Operations</b>	<b>12</b>
3.1	Introduction	12
3.2	First use	12
3.3	Important instructions	12
3.4	Operating the fitvibe smart	12
3.4.1	Operating the apparatus in general	12
3.4.2	System settings	13
3.4.3	Summary of operating modes	13
3.4.4	Start-stop switch	14
3.4.5	Use of hand straps.	14
3.4.6	Vibration training parameters	14
3.4.7	Increase training effort	15
<b>4</b>	<b>Technical problems, warranty, maintenance &amp; cleaning</b>	<b>16</b>
4.1	Technical problems	16
4.1.1	The device does not want to start.	16
4.1.2	No function seems to be working	16
4.2	Warranty and service	16
4.3	Maintenance	16
4.4	How to clean the device	16
<b>5</b>	<b>Technical Data</b>	<b>17</b>

# 1 Safety

## 1.1 Purpose

The fitvibe smart is intended solely for vibration training. Before attempting to use the fitvibe smart, read and understand the safety aspects and user instructions in this manual.

## 1.2 Legal aspects of training

Training on the fitvibe smart is intended only for those who are physically capable of undertaking a course of exercise. The legal regulations of the country apply.

## 1.3 Safety instructions

### 1.3.1 General

- Only a technician authorized by Uniphy Elektromedizin GmbH & Co. KG is entitled to carry out repairs or modifications and to add elements to the device.
- The electrical system of the user must comply with all laws and regulations.
- Follow the instructions and the directions in these user instructions.
- Always keep these user instructions with the equipment.
- Place the equipment on a horizontal and stable base.
- Do not place any objects on the equipment.
- Do not place the equipment in the sun or above a heat source.
- Only use the device at a temperature in between +10°C and +40°C.
- Do not use the equipment in a humid environment.
- Do not let any liquid flow into the equipment.
- Do not disinfect or sterilize the equipment; Clean the equipment with a dry or moistened cloth.
- All legal restrictions concerning the use of the device must be complied with.

### 1.3.2 Contra-Indications

Each user must check whether he or she suffers or might be suffering from one or several of the following contra-indication before starting the vibration plate training. If so only use the fitvibe smart when you have approval from a medical specialist.

<b>Absolute contraindications:</b>	<b>Relative contraindications:</b> (use fitvibe only after consulting a medical doctor.)
Acute inflammations, infections and/ or fever Acute arthropathy and arthrosis Acute rheumatoid arthritis Acute headache Postoperative wounds Implants / Endoprotheses Metal or synthetic implants (e.g.. pacemaker, artificial cardiac valves, recent stents) Acute thromboses or increased thrombotic afflictions (6 months) Acute back problems Pregnancy Severe osteoporosis with BMD < 70mg/ml Spasticity Morbus Sudeck Stadium I (CRPS I) Tumors with metastases in the musculoskeletal system Benign Paroxysmal Positional Vertigo (BPPV)	Epilepsy Gallstones, kidney stones, bladder stones Articular rheumatism and arthrosis Heart failure, insufficiency NYHA III-IV Cardiac dysrhythmia Cardiac affections Metal or plastic implants (piercing or gynecologic implants Chronic back pain (6 months after acute appearing) Severe diabetic mellitus with arterial circulation disorders Tumors (in a period 5 years plus after appearing, excluding metastases in the musculoskeletal system) Spondylolisthesis without gliding M. Parkinson Chondromalasy of the joints of the lower extremities, osteonecrosis and -chondrosis Arterial circulation disorders (Stadium Fontaine III°) Venous insufficiency with ulcus cruris Morbus Sudeck Stadium II (CRPS) Lymphatic edema

### 1.3.3 Unexpected symptoms

A user must immediately stop using the fitvibe smart when any unexpected symptoms occur. Examples are:

- Muscle contractions without control during static exercises.
- Irregular movement of the device during dynamic exercises.
- Dizziness
- Pain
- Muscle tightness
- Nausea

### 1.3.4 Electrical Safety

- Only use the equipment in an area with facilities that meet the applicable legal regulations.
- Connect the equipment to an outlet with a protective earth terminal. The outlet must meet the locally applicable requirements.
- Fitvibe smart has been manufactured for use at a power voltage of either 230V, 50-60Hz or 115V, 50-60Hz. Check the type shield for compatibility. Use of the device at an unsuitable voltage can cause irreparable damage

### 1.3.5 Prevention of explosion

- Do not use the equipment in an area where combustible gases or vapours are present.
- Switch off the equipment when it is not used.

### 1.3.6 Liability

Neither the manufacturer nor the local dealer can be held liable for personal injuries or for damage to or by the equipment used, if for example;

- The equipment or the accessories are used incorrectly
- The user instructions are wrongly interpreted or ignored
- The equipment is badly maintained
- Maintenance or repairs are performed by people or organisations that are not authorized to do so by Uniphy Elektromedizin.

### 1.3.7 Conformity

MDD/CE

The fitvibe smart has been developed and is tested, produced and sold according to the international quality standards of the ISO 9001:2000 and DIN EN ISO 13485:2003 quality assurance system. Furthermore, the fitvibe smart meets the strictest requirements of European Directive 93/42/ECC for Medical Devices (MDD), as indicated by the presence of a CE identification (CE 0197) on the apparatus. Furthermore, the certifying authority (notified body), TÜV Rheinland/Berlin-Brandenburg, also carries out checks regularly. This lets you, the treating physician, be certain that you are working with equipment that meets all the requirements and that provides you and your patients with the best possible levels of safety. For the manufacturer's liabilities, see Chapter 4.

### 1.3.8 Deactivation and disposal

Observe the following instructions when deactivating the "fitvibe smart" vibration device:

You may deliver the device including all accessory items to a waste disposal service provider specialized on electronic products or capable of forwarding the device to such a specialist. In countries where the suppliers are under the obligation of taking back deactivated electronic devices, you may of course contact your supplier. Most countries have by now introduced legislation which prevents consumers from simply throwing their deactivated electronic products away alongside their ordinary household waste. In the event that you yourself are responsible for the safe disposal of the system or that you prefer to do it on your own, we would like to provide you with the following information:

The components of the "fitvibe smart" vibration device can be classified as follows: Vibration plate, handle with electronic unit, cable, straps and soft mat must be treated like electronic scrap or light chemical waste. These components contain tin, copper, iron, a range of other metals and different plastics. Most countries regard these types of products as light chemical waste. Please inform yourself about the waste disposal regulations which apply in your country. Packaging materials and instructions can be recycled. You may want to bring them to an appropriate collecting point or add them to your household waste, depending on how your local waste disposal system has been organized.

Everything said in this paragraph about the safe disposal of electronic articles applies equally to the accessory items.

## 2 Components and Installation

### 2.1 *Receipt*

1. Check whether the equipment has been damaged during transport.
2. Check whether the accessories are intact and complete.
  - Inform your supplier immediately of any damage, defects or incomplete delivery.
  - Do not use the equipment in case of damage or defect!

## 2.2 Description of the control unit and accessory items

In Figure 1 you see a diagram of the fitvibe smart with all accessory items. The numbers in the diagram and in the list below refer to the same items.

1. Vibration plate
2. Hand straps (2 pieces)
4. Handle
6. Control Unit/ display
7. "Start/Stop" button
8. Button "Arrow UP"
9. Button "OK"
10. Button "Arrow DOWN"
11. Display
16. Side cushion (optional accessory)
17. Soft Mat (1 piece)
18. Power cable
19. Power switch
20. Fitting for straps
21. Socket to connect remote control
22. Remote control "Start/Stop"

Besides that in every package there is a

- Poster with exercises
- Fitvibe for dummies
- Manual



Figure 1

## 2.3 Assemblage

1. Place the fitvibe smart on a horizontal and stable base.
2. Insert the handle (4) in the holes of the control unit and fix with 2 screws at each side.

3. Use the designated connection (21) to connect the remote control (22).
4. Check whether the mains voltage is that is stated on the type shield corresponds with the voltage of your mains supply (see 1.3.4)
5. Connect the power connection (18) of the device to an outlet
6. Switch the power switch (19) on.

## 2.4 Transport & Storage

Take account of the following matters if the equipment has to be transported or stored:

- Transport or store the equipment in the original packaging
- The maximum period for transport or storage is 15 weeks
- Temperature: -20°C to +60°C
- Relative humidity: 10% to 100%
- Atmospheric pressure: 500hPa to 1060hPA

## 3 Operations

### 3.1 Introduction

The "fitvibe smart" generates three-dimensional mechanical vibrations of the main vector in vertical direction.

The "fitvibe smart" is extremely easy to operate. The operating principle is clear and self-evident.

### 3.2 First use

Make sure that a qualified coach supervises the first training session.

### 3.3 Important instructions

- Check that no contra-indications are present. (see pt 1.3.2).
- Stop using the fitvibe smart immediately when any unexpected symptoms occur. (see pt 1.3.3).
- Wear shoes with a flexible sole during training.
- Use the soft mat when you support your weight other than with your feet.
- Follow the instructions for the exercise and the rest time to prevent injuries.
- Start by standing on the tiptoes with the knees unlocked to absorb the maximal amount of the vibration through your calves and thighs.
- Do not extend the limbs that are in contact with the vibration plate. Correct your position if you feel vibration in your head, voice and/or eyes being are the best indicator of a wrong position.
- Take a sports drink or water, approximately 30 to 60 minutes, before the start of the training session.

### 3.4 Operating the fitvibe smart

#### 3.4.1 Operating the apparatus in general

- By means of the UP and DOWN arrow (8,10), one can scroll trough the menu.
- Press the OK button (10) to activate the displayed menu item
- Change settings via the UP and DOWN arrow (8,10)
- Press the OK button (10) longer than 0,3 sec to return one step in the menu.


Figure 3

- When the power is switched on, a welcome screen is displayed giving 3 possibilities:
  - Press the arrow UP button (8) for more than 0,3 sec and the list of contra-indications will be displayed (see pt.1.3.2)
  - Press the arrow DOWN button (10) for more than 0,3 sec and the system settings will be displayed (see pt 3.4.2)
  - OK button (9) continue and indications menu is displayed (see pt 3.4.3)

### 3.4.2 System settings

The menu items are:

- **Language**  
Text display language: English, Dutch, German, ...
- **Volume of "beep"**  
Modifiable between 0% (no sound) and 100% (maximal sound).
- **Button beeps**  
Yes: Sound as buttons are pressed is switched on.  
No: Sound as buttons are pressed is switched off.
- **Time up beeps**  
Yes: Audio alerts announce the end of an activity or rest period.  
No: No audio alerts are generated to announce the end of an activity or rest period.
- **Active time in the Free quick start programme**  
The standard given time in the quick start menu can be modified in different steps between 10-180 sec.
- **Intensity in the Free quick start programme**  
The standard given vibration intensity van be modified in 6 levels, influencing amplitude and frequency of vibration)

### 3.4.3 Summary of operating modes

The main menu enables the selection of several indications and a free mode.

- **Indications**

Indications are default programs in the fitvibe smart. For every indication there are 4 levels of training schedules with different exercises and parameters.

Pictures show clearly how to execute the correct position, clicking the OK button (9) shows the training parameters. (The display also switches back automatically after a couple of seconds).

You can browse forward and backward through the exercises by using the Up and Down arrows

The word 'Exercise' is displayed on the screen's header. The exercise's number within the sequence is shown to the right.

Example: Exercise 3 / 5: The exercise currently running is the third of a total of five. The exercise number is displayed top left below the header, e.g. 37.

While the active time is running down and during any pause, a different time screen is shown.

Whenever active time is running, both this and the associated pause time can be deleted with a very long press of the STOP button.

If no active time or pause is running, a very long press of the OK button will take you from the time screen back to the parameters menu and via another long press on the OK button you return to the main menu.

#### ▪ **Free Modus**

In this menu one can start the vibrations immediately with the preset parameters by clicking the start/ stop button (7). With the button arrow UP (8) and DOWN (10) the user is also able to change time and training intensity.

While the active time is running down, a time screen with extra large figures is displayed (active time display).

Whilst training training intensity can be changed with the button arrow DOWN (10)

When the time is stopped (with the start/ stop button) one can always return to the parameter menu via a long press on the OK button.

### 3.4.4 Start-stop switch

The "fitvibe smart" vibration device has been equipped with two start-stop switches. Both are capable of starting and stopping the vibration independently from one another:

One orange-coloured start-stop button (7) are located on the control unit; the remote control (22) is connected by cable with the control unit.

### 3.4.5 Use of hand straps.

Straps can be used to transfer the vibration to certain sections of the body as shown in the illustration. The non-elastic straps serve to convey the vibrations to the muscles in question. The straps can be fastened to the vibration plate at the fitting for straps (20). The straps have been provided with reference lines in different colours, allowing the user to adjust the length of one strap to the length of the other.

### 3.4.6 Vibration training parameters

Vibration frequency:	30 – 35 - 40 Hz		
Intensity (amplitude):	high / low		
Active time:	0 - 180 sec		
Operating modes:	Relax & Massage		(4 performance levels)
	Upper body training	(" "	")
	Lower body training	(" "	")
	Stretching	(" "	")
	Body Shaping	(" "	")
	Free Mode		

### **3.4.7 Increase training effort**

You may increase the required effort in a number of different ways such as:

- Increase the vibration frequency
- Extend the active periods
- Shorten the rest periods
- Increase the intensity from low to high
- Increase the number of repetitions

You can also increase the required effort by making the exercises more difficult to perform. You can achieve this by:

- selecting more difficult / more complicated starting positions (see poster)
- change from static to dynamic exercises
- enlarge the stretched muscle section

## **4 Technical problems, warranty, maintenance & cleaning**

### **4.1 Technical problems**

#### **4.1.1 The device does not want to start.**

Check whether there is any power, whether the power switch is on and whether power cable and fuses are okay. If the fuses have blown, you may want to replace them yourself, but make sure that you only use the type of fuse indicated on the type plate. You can ask your supplier to order the right fuses for you.

#### **4.1.2 No function seems to be working**

Switch the device off and on again. Contact your supplier if the problem persists.

### **4.2 Warranty and service**

- Only duly authorized service engineers are entitled to carry out repairs or modifications and to add elements to the device.
- The electrical system of the user must comply with all laws and regulations.
- The device must only be operated by duly authorized persons in strict compliance with the user manual.
- The device must only be used for its designated purpose.
- The device must be regularly maintained in compliance with the maintenance instructions.
- All legal restrictions concerning the use of the device must be complied with.

Apart from the fuses, there are no components inside the "fitvibe smart" device which may be exchanged and replaced by the user. In the event of incorrect use and non-compliance with the maintenance instructions, the Uniphy Elektromedizin GmbH & Co. KG and their duly appointed representatives will cease to be liable for all resulting damages, personal injuries, defects or malfunctions. Customer service and warranty will be provided through your supplier. The terms and conditions of your supplier will apply. The warranty will be invalidated by any violation of the aforementioned instructions. The device is expected to have a technical life of 10 years.

Fitvibe SMART is mentioned for a "non professional environment"

Guaranteed period: 2 years with a maximum of 1000 operating hours. The average use of the fitvibe SMART should not exceed 1 hour (vibration time) per day. A more intensive use is considered by the manufacturer to be consistent with a professional environment for which the fitvibe Medical and Fitvibe Excel series have been developed.

### **4.3 Maintenance**

The manufacturer recommends subjecting the device to an annual inspection by an engineer. The following components should be checked with particular attention:

- all bolted and mechanically clamped connections
- all cables and their connecting elements (visual control)
- the functional safety of the manual switch and the foot switch
- the rubber springs (visual check for mechanical damage)

### **4.4 How to clean the device**

You can use a slightly moist cloth to clean the vibration plate, the control unit and the accessory items. Use a mild and non-abrasive detergent!

The handle, if necessary, can be treated with a commercial disinfectant.

## 5 Technical Data

The "fitvibe smart" is a Whole Body Vibration (WBV) device.

Vibration frequency:	30 - 35 - 40 Hz
Intensity (Amplitude):	high / low
Active time:	0 - 180 sec
Power voltage:	115 Volt $\pm$ 10 %, 50Hz/60Hz or 230 Volt $\pm$ 10 %, 50Hz/60Hz
Dimensions (W x D x H):	630 x 940 x 1455 mm
Weight:	ca. 70 kg
Ceiling load:	200 kg max. static load and a dynamic proportion of 20 %
Noise level:	max. 69 dBA

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