THERMO 500
SHORTWAVE THERAPY UNIT

USER MANUAL
GEBRUIKERSHANDLEIDING
GEBRAUCHSANWEISUNG
MODE D'EMPLOI

EN - NL - DU - FR

Your supplier is:

CE 0197

--

Uniphy Elektromedizin GmbH & Co. KG
Neuendorfstrasse 19 b
D-16761 Hennigsdorf - Germany
T +49 (0) 3302 5044-0
F +49 (0) 3302 5044-99

World wide distributor:
GymnaUniphy NV
Pasweg 6A
3740 Bilzen
Belgium
T +32 (0) 89 510.510
F +32 (0) 89 510.511
info@gymna.com
www.gymna.com

ENGLISH.....	11
NEDERLANDS	43
DEUTSCH.....	77
FRANÇAIS	109

INTRODUCTION.....	11
1 OPERATION ELEMENTS AND SYMBOLS	12
1.1 INTRODUCTION	12
1.2 OPERATION ELEMENTS AND SYMBOLS	12
2 IN GENERAL.....	15
2.1 USE OF HIGH-FREQUENCY HEAT THERAPY	15
2.2 THE APPLIANCE	15
2.3 ELECTRODE ARMS	15
2.4 ELECTRODES.....	16
2.5 CONTACT-CONTROL IN CASE THE DISTANCE BETWEEN THE THERMOPLODES AND THE SKIN IS TOO BIG.....	16
2.6 EUROPEAN DIRECTIVE FOR MEDICAL DEVICES 93/42/EEC - MEDICAL DEVICES DIRECTIVE (MDD)	17
3 OPERATION	18
3.1 CHECK ON DELIVERY	18
3.2 CURRENT CONNECTION	18
3.3 SWITCHING ON.....	18
3.4 FUNCTIONAL TEST	19
3.5 SELECTION OF LANGUAGE	19
3.6 SENDING / SHIPPING (TRANSPORT).....	19
4 SELECTION OF THE SETTING MODE	20
4.1 SELECTION VIA THE STANDARD MODE	20
4.1.1 <i>There is no Thermoplude connection</i>	20
4.1.2 <i>There is one Thermoplude connection</i>	20
4.1.3 <i>Two Thermoplodes 140 are connected</i>	21
4.1.4 <i>Two different Thermoplodes are connected</i>	23
4.2 SELECTION VIA THE INDICATION LIST.....	25
4.2.1 <i>End of an indication programme</i>	26
4.3 SELECTION VIA THE FREE MODE	28
4.3.1 <i>Characteristics of the cursor</i>	29
4.3.2 <i>Creation of a new programme</i>	29
4.3.3 <i>Execution of the programme</i>	34
5 CONTROL OF THE ELECTRODE CONNECTION.....	35
6 TREATMENT	36
6.1 PREPARING THE PATIENTS	36
6.2 THERMAL EFFECT OF THE COIL-FIELD METHOD.....	36
6.3 INDICATIONS.....	36
6.3.1 <i>Contra-indications</i>	37
7 MALFUNCTIONS, GUARANTEE, RELIABILITY, CUSTOMER CARE	38
7.1 MALFUNCTIONS.....	38
7.2 GUARANTEE AND RELIABILITY	38
7.2.1 <i>Terms of guarantee</i>	38

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INTRODUCTION

Congratulations!

You have opted for the Thermo 500, a short-wave therapy appliance with splendid performances and a nice design, general-purpose practical and with a permanent user's quality.

In physiotherapy practices, at the doctor's and in hospitals the Thermo 500 has already proved its achievements.

In the development of this appliance we have attached great interest to a high degree of reliability, the safety, user-friendliness and a long life span.

To use the Thermo 500 correctly, please read the instruction manual first.

We wish you and your patients every success with the treatment using the Thermo 500.

GymnaUniphy NV

1 OPERATION ELEMENTS AND SYMBOLS

1.1 Introduction

On the following page of this manual, you can find a drawing of the appliance. The numbers mentioned below correspond to the numbers on the drawing.

1.2 Operation elements and symbols

- | | | |
|--------|--|---|
| 1. | Mains switch | |
| 2. | Start/stop key |  |
| 3. | Standard mode key 'Standard' |  |
| 4. | Free mode key 'Free mode' |  |
| 5. | Indication mode key 'Indication' |  |
| 6. | 'Select <' key |  |
| 7. | 'Set/power' turning knob | |
| 8. | LED 'set', setting parameters active | |
| 9. | LED 'power', setting power active | |
| 10. | 'Select >' key |  |
| 11. | LED 'Select active' | |
| 12. | Display | |
| 13. | Cable electrode fixing point | |
| 14. | Adjustment regulator to fix the electrodes | |
| 15. | Hinge | |
| 16/17. | Locking/Fixing knob | |



- 18. Bottom part of the electrode arm
- 19. Locking/Fixing knob to turn loose/ fix the vertical movement of the arm
- 20. Fastening bracket
- 21. Bottom rotation hinge of the arm
- 22. Screws to fix the electrode arm
- 23. Contact points for electrode cable
- 24. Mains connection / Fuses


Please read the instruction manual first.

Not ionized radiation

Electrode arm

Back side of the device

2 IN GENERAL

2.1 Use of high-frequency heat therapy

The Thermo 500 is a high-frequency electrotherapy appliance, functioning on the set frequency of 27,12 MHz (wavelength 11m). It enables using high-frequency heat therapy in the coil-field, as well with one channel as with two channels. The appliance is operated in the pulse mode, that is why it is suitable for all thermal and a-thermal applications. The use of high-frequency energy for heat therapy offers the advantage of a larger infiltration depth, compared to simpler procedures such as compresses, baths, infrared light, thermal pillows and microwaves.

Endogenous heat causes a number of physiological processes affecting muscles, thews and other interstitial tissue structures in a paroxysm removing manner, it increases the cellular metabolism, the enzyme quickness of reaction and a continuous blood circulation in the treated part.

With the possibility to apply the high-frequency energy in short and high energy pulses (pulsed working) the depth working, in particular the increasing effect on the blood circulation, can be raised further while a heat stimulus can hardly be felt on the heat-sensitive skin.

The application scope of high-frequency heat therapy is very extensive. Into prominence there are all rheumatic disorders and pains of muscles and joints, inflammation diseases of the respiratory organs, the kidneys and urinary passages and all diseases caused by a bad blood circulation.

2.2 The appliance

The Thermo 500 is an easily rolling device on large swivel castors. By means of the break on the two rear castors you can prevent it from rolling when it is not desirable. On the left of the appliance the mains switch **(1)** is installed. Setting the power is done via the turning knob **(7)**.

2.3 Electrode arms

At the reverse of the appliance, amongst others, the screw holes are situated to fix the electrode arm(s) **(22)** and the connection contacts **(23)** for the cable of the coil-field electrode(s), which are called "Thermoplodes" for the Thermo 500. The electrode arms are assembled at the reverse of the appliance by means of the brackets **(20)**.

Through its stable construction and five hinges, the electrode arm will be able to fix the Thermoplodes steadily in all kinds of treatment positions. Except for the horizontal turn of the arm and the hinge **(15)** in the Thermoplode fixing point **(13)** the user can modify all the hinges according to specific needs. The fixing knob **(19)** of the bottom rotation hinge **(21)** makes it possible to secure the entire electrode arm with the arm maximally extended and with different electrode weights.

After turning the locking knobs loose by means of the other main hinges **(16)** and **(17)** length and height variations can be applied in a mutually combined action. Herein the electrode arm needs to be supported at the upper hinges **(15)** and **(16)**.

2.4 Electrodes

The Thermo 500 is an appliance for one- or two-channel short-wave therapy. Standard, 1 coil-field electrode with a 14 cm diameter is supplied with the appliance (see also chapter 10 *SURVEY OF THE ACCESSORIES*).

The Thermoplodes and cables are especially adapted to the Thermo 500. So please do not use other cables and electrodes for safety reasons.

2.5 Contact-control in case the distance between the Thermoplodes and the skin is too big

When the contact with a Thermoplude becomes insufficient as a result of a movement of the patient, (if the electrode is not placed on the skin)) the display shows this. Nevertheless, the treatment time continues. In case the situation is still the same after 10 seconds, there is an additional whistle alarm and the treatment time stops. After an adjustment of the contact, the treatment will automatically be continued.

Caution!

1. For an optimum effect, the Thermoplude has to be placed directly on the area to be treated.
2. To set the power please pay attention to the instructions for use of the Thermo 500 and the parameter proposal of the indications. In particular you should pay attention in case of Thermal treatments. The feeling of warmth will have a delaying effect on structures which are deeper.
3. You should adjust the Thermoplude to the size of the area to be treated. The surface to be treated should correspond to the surface area of the electrodes. If the contact is bad or in bony areas the message “improve position Thermoplude” appears.
4. If there is no contact and the power is set to high, the Thermoplude will heat up internally, in particular near to the cable connection.
5. Defects of the Thermoplude due to incorrect use are not covered by the warranty.

When cleaning and disinfecting the appliance, prevent any liquid from entering the appliance or parts of the accessories. Contact points, which have become wet, need to be dried thoroughly before further use.

2.6 European Directive for Medical Devices 93/42/EEC - Medical Devices Directive (MDD)

The Thermo 500 meets the MDD- requirements. It complies with the general demands of labour safety by applying the valid technical standards and observing the conditions in technical files.

According to the MDD this electro-medical device can only be used by persons who can guarantee a professional application based on their education or their knowledge and practical experience and who have been trained in order to use the appliance correctly taking into account the instructions of this manual.

As manufacturer we can only be responsible and liable for the technical safety characteristics of the appliance, when the Thermo 500 is used in accordance with the instructions of the manual. Repairs, also opening the appliance, can only be executed by us or by customer service partners authorised by us.

For the sake of the patients' and users' safety we recommend an annual safety test. We have already experienced that most of the assumed malfunctions can often be linked to mistakes and faults in using the appliance. So please check on the operation of the appliance before calling the customer service.

Whoever using the appliance, should check on the operation and the condition thereof. Regular check-ups of all the cables and wires on isolation defects are essential.

Caution!

The appliance arouses high-frequency electric and magnetic fields coming right through walls, ceilings and floors. It is inevitable that part of these fields occur near the appliance.

Electronic devices in the vicinity of the Thermo 500 can be affected. The disturbance strongly depends on the distance between the devices. Please notice to place the appliance at a distance of at least 5 metres from surrounding electronic (interference-sensitive) equipment and pay attention not to point the Thermoplodes in the direction of other electronic devices such as electro-stimulation equipment and other electronic apparatus. This problem can be solved completely when the short-wave appliance is put in a protected area, meaning a Faraday cage (a Faraday cage does not allow electromagnetic fields).

3 OPERATION

Caution!

The appliance cannot be used in an environment endangered with explosions. When explosive narcotics are used simultaneously in the anaesthetic room, the danger of an explosion cannot be excluded.

3.1 Check on delivery

Check whether the appliance has not been damaged during transport and whether the accessories are intact and complete (*see chapter 10 SURVEY OF THE ACCESSORIES*).

In the event of damage or defects, contact your supplier. Do not switch the appliance on in case of serious damage. When the appliance is damaged after the first use, you should have it checked by an authorised organisation.

3.2 Current connection

The appliance must be connected to a mains voltage of $230\text{ V} \pm 10\%$ and 50/60 Hz. There is a special edition for $115\text{ V} \pm 10\%$ and 50/60 Hz. Mind the model plate on the reverse of the appliance. Before connecting to the current network, you must check the accordance with the model plate.

The connection can only happen by means of the current cable being part of the basic equipment and a regulatory installed socket. For the fuse we recommend 10 A with 230 V.

3.3 Switching on

You can switch the appliance on by means of the mains switch **(1)** on the left-hand side of the console. The name of the appliance and the software version are shown for a while on the display. During an internal test of the appliance all the blinkers (indicating lights) start to illuminate shortly and after a buzzer sound the appliance is set in the standard mode (the standard mode LED lights up and the concerned information is shown on the display).

					T	h	e	r	m	o		5	0	0						
				V	e	r	s	i	o	n		1	1	.	x	x		x		

The following settings are done via the console.

3.4 Functional test

The appliance is tested during production for electrical safety. Each time the appliance is switched on, the microprocessor carries out a profound functional test to assess afterwards whether the displays function correctly. If a fault is detected, you cannot use the appliance and you must have it repaired.

3.5 Selection of language

The appliance offers various language options. When you keep the “SELECT ◀” and “SELECT ▶” - keys pressed while switching on, the Thermo 500 enters the setting mode for the user’s language.

						E	N	G	L	I	S	H								

The language can be selected by means of the turning knob: English, Nederlands, Deutsch or Français.

By pressing the start/stop key  the appliance is ready-for-use.

3.6 Sending / shipping (Transport)

In case the appliance needs to be sent/transported, we advise you to use the original packaging.

4 SELECTION OF THE SETTING MODE

4.1 Selection via the standard mode

When switched on the Thermo 500 will automatically start in the standard mode. As the various parameters have already been pre-programmed, you can immediately start the treatment in this mode. You only have to set the desired power and possibly the treatment time.

Also after pressing the "standard mode key"  the Thermo 500 comes into standard mode. The LED near the key lights up.

The parameter values depend on the connected and automatically traced Thermopode(s). The various situations are described hereinafter.

4.1.1 There is no Thermopode connection

						S	T	A	N	D	A	R	D					
								2	0	m	i	n						

The programme cannot be started. After the start/stop key  is pressed the Thermo 500 alarm buzzer will sound.

Immediately after the Thermopode is connected the concerned electrode symbol appears on the display and the programme can be started.

4.1.2 There is one Thermopode connection

A Thermopode 140 is for example connected to the left channel, the START window looks as follows:

						S	T	A	N	D	A	R	D					
7	0	.	0	W				2	0	m	i	n						

In the START window you can see the connected Thermopode, the average power and the treatment time.

The parameters of the connected Thermoplode 140 are :

pulse time 400 μ s
pulse frequency 875 Hz
average power 70 Watt

The parameters of a connected Thermoplode 80 are:

pulse time 400 μ s
pulse frequency 500 Hz
average power 32 Watt

The parameters if a Thermoplode 80 and a Thermoplode 140 are connected simultaneously are:

Thermoplode 80	Thermoplode 140
pulse time 400 μs	pulse time 400 μs
pulse frequency 500 Hz	pulse frequency 500 Hz
average power 32 Watt	average power 40 Watt

4.1.3 Two Thermoplodes 140 are connected

As well to the left as to the right channel a Thermoplode 140 is connected, the START window looks as follows:

						S	T	A	N	D	A	R	D						
4	0	.	0	W				2	0	m	i	n		4	0	.	0	W	

4	0	0	μ	s			5	0	0	H	z				4	0	0	μ	s	
	2	0	0	W	□			2	0		m	i	n			2	0	0	W	□

The previously set treatment time is 20 minutes. This is the only parameter which can be set now by means of the turning knob (the setting is also possible in the START window by means of the standard mode). The parameter flashes the moment the treatment time can be set.

					S	T	A	N	D	A	R	D							
0	.	0	0	W				2	0	:	0	0		0	.	0	0	W	

The start value of the average power is 0.00 W (and the peak power 0 W). The cursor flashes when the power can be set. The average power of the left channel can be increased by the turning knob. After starting the treatment, the treatment time parameter counts back in steps of 1 second.

					S	T	A	N	D	A	R	D							
4	0	.	0	W				1	9	:	5	0		0	.	0	0	W	

Press the SELECT > key to select the right-hand channel. The indication of the average power of the right channel flashes and can now be increased by means of the turning knob.

					S	T	A	N	D	A	R	D							
4	0	.	0	W				1	9	:	4	0		4	0	.	0	W	

Press the standard mode key to display the parameters. The parameter settings appear on the display. This window/screen displays the set peak power value.

4	0	0	μ	s				5	0	0	H	z				4	0	0	μ	s
2	0	0	W					1	9	:	4	0				2	0	0	W	

At the end of the treatment time, you can hear a buzzing sound of about 10 seconds and the position of before the start is recovered, meaning set back to the initial value; the power is automatically set back to 0 W.

Treatment can also be stopped before the end of the treatment time by pressing the

start/stop key .

4.1.4 Two different Thermoplodes are connected

A Thermoplode 140 has been connected to the left-hand channel and a Thermoplode 80 to the right-hand channel. The START display looks like this:

						S	T	A	N	D	A	R	D						
4	0	.	0	W				2	0	m	i	n		3	2	.	0	W	

To display the parameters, press the standard mode key again

If a Thermoplode 80 is connected you should bear in mind that the average power may not exceed 32 Watt. The Thermoplode 140, which is connected to the other channel, then has a maximum average power of 40 Watt.

4	0	0	μ	s				5	0	0	H	z			4	0	0	μ	s
2	0	0	W					2	0	m	i	n			1	6	0	W	

The preset treatment time is 20 minutes. This parameter is the only one, which can be set at this time with the turning knob (it can also be set in the START window of the standard mode). The parameter flashes to indicate that the treatment time can be set.

Treatment can be started with the aid of the start/stop key . The START window is again shown on the display :

						S	T	A	N	D	A	R	D						
0	.	0	0	W				2	0	:	0	0		0	.	0	0	W	

The starting value of the average power is 0.00 W (and the peak power 0 W). The cursor flashes if the power can be set. The average power of the left-hand channel can be increased with the turning knob. When treatment starts the treatment time will count down in steps of 1 second.

						S	T	A	N	D	A	R	D						
4	0	.	0	W				1	9	:	5	0		0	.	0	0	W	

Press the SELECT >key to select the right-hand channel. The indication of the average power of the right-hand channel flashes and can at this time be set with the aid of the turning knob.

					S	T	A	N	D	A	R	D							
4	0	.	0	W				1	9	:	4	0		3	2	.	0	W	

Press the standard mode key  to display the parameters. The parameter settings will appear on the display. The peak power set is shown in this screen.

4	0	0	μ	s				5	0	0	H	z				4	0	0	μ	s
2	0	0	W					1	9	:	4	0				1	6	0	W	

When you press the start/stop key , the treatment starts. At the end of the treatment time, you can hear a buzzing sound of about 10 seconds and the position of before the start is recovered, meaning set back to the initial value; the power is automatically set back to 0 W

4.2 Selection via the indication list

The Thermo 500 disposes of a very extensive list of pre-programmed indications with advice on the treatment parameters. In this way you can define your setting quickly and simply, as well for well-known as for unknown pathologies. By means of the centrally placed turning knob you can easily page the indication menu and make your choice. The indication menu is drawn up in a way that takes into account the current situation of the disorder. In the indication mode only the parameters treatment time and average power can be set (up to the maximum peak power).

Caution!

For indications, the Thermo 500 automatically detects the type of Thermoplane that is being connected. You can choose which type of Thermoplane you wish and to which channel you want to connect the Thermoplane.

If you opt for a indication which has a higher average power than 32 Watt and you have connected a Thermoplane 80, the Thermo 500 will automatically reduce the average power to 32 Watt. This will reduce the peak power.

In that case a flashing Thermoplane 140/80 symbol will indicate that you have connected the Thermoplane 80 (see chapter 5).

By pressing the indication mode key  the Thermo 500 comes in **indication mode**. The LED of the indication key lights up. The last used indication appears on the display; after switching on the appliance, you always see the first indication of the list. The list is arranged in alphabetical order.

	C	L	A	U	D	I	C	A	T.		I	N	T.				
6	.	0	0	W				2	2	m	i	n					

The appliance contains 38 indications (available in 4 languages). By means of the turning knob you can page the indication list.

	A	R	T	H	E	R	I	T	I	S							
2	.	4	7	W				1	6	m	i	n					

To see the parameters you have to press, from the START window illustrated above, the indication mode key  again. The parameter settings appear on the display:

2	0	0	μ	s			6	5	H	z						
	1	9	0	W			1	6	m	i	n					

To return to the START window, you press the indication mode key  again.

4.2.1 End of an indication programme

In the indication mode only the parameters treatment time and average power can be set. After you have selected the indication, you will see the START window with the advised parameters on the screen (the advised average power, in the example below 2.47 W, is the maximum value which you can set after pressing the start/stop key).

Press the SELECT > key  to adapt the treatment time. The value flashes; by means of the turning knob you can change the treatment time. (this can also be done in the parameter window):

	A	R	T	H	E	R	I	T	I	S						
2	.	4	7	W			1	2	m	i	n					

By pressing the start/stop key , the treatment is started at the initial value of 0 W of the average power.

	A	R	T	H	E	R	I	T	I	S						
0	.	0	0	W			1	2	:	0	0					

By means of the turning knob the average power of the left channel can be raised to the maximum average power of the displayed programme. The time indication counts back in steps of 1 second.

	A	R	T	H	E	R	I	T	I	S						
2	.	2	0	W			1	1	:	5	0					

When the electrode on the right-hand side is also connected, you can raise the average power of the right channel by pressing the SELECT > key  to the maximum peak power of the displayed programme.

If you wish to see more information on the parameters, you have to press the

indication mode key . Subsequently the parameter settings appear on the display. The average power can also be set in this window.

At the end of the treatment the treatment is terminated and the parameters are set back to the initial values of the programme. At the end of the treatment time you will

also hear a buzzing sound of about 10 seconds. When you press the  button during treatment the power will go to “0” and the treatment time will give the time when stopped in minutes

4.3 Selection via the free mode

When you press the free mode key  the Thermo 500 passes to **free mode**. On the display “FREE DESIGN” appears and the LED near the key lights up. The free mode contains 50 programme places, which can be freely programmed and saved, as you wish. Only the first programme is a FREE DESIGN programme, which cannot be saved.

The Thermo 500 automatically traces which Thermoplode is connected.

				F	R	E	E		D	E	S	I	G	N					
								2	0	m	i	n							

By means of the turning knob you can page the programmes. The pre-programmed names of the programmes are FREE 01 to FREE 49.

In the following example the Thermoplode 140 is connected to the left channel in the free mode FREE 03. The previously set treatment time is 20 minutes.

					F	R	E	E		0	3							
0	.	0	0	W				2	0	m	i	n						

When you wish to see the parameters of this programme, you have to press the free mode key  again. Subsequently the parameter settings appear on the display.

The following parameter values are standard pre-set for all programmes.

Pulse time = 65 μ s
Pulse frequency = 25 Hz
Peak power = 10 W
Treatment time = 20 minutes

	6	5	μ	s				2	5	H	z					6	5	μ	s
	1	0	W					2	0	m	i	n				1	0	W	

Parameters changed by you are automatically saved, even when the treatment has not been started yet by means of the start/stop key .

If you change parameters after pressing the start/stop key  and after starting the treatment, these changes (parameters adapted during the treatment) are not saved. After the treatment time or after stopping the treatment by pressing the start/stop key  the initial values appear again.

4.3.1 Characteristics of the cursor

You can change the position of the cursor on the display (the flashing character which can be set by the turning knob), by pressing the SELECT  or SELECT  key. Whenever a field is empty, an underlining flashes. By pressing SELECT the fields are addressed in the following order.

4.3.2 Creation of a new programme

Overwriting a free programme or a previously saved programme creates new programmes.

Basic principle of programming:

- Select the parameter to be set by pressing SELECT  or SELECT 
- Set the desired value of the selected parameter (flashing now) by means of the turning knob
- Press SELECT  or SELECT  to pick out the following parameter

Caution!

When a parameter is changed, it is saved automatically. Pay attention not to overwrite an existing programme unintentionally.

From the moment a programme can be started (this is when all relevant parameters of at least one channel are set correctly), the LED key next to the start/stop key  lights up.

4.3.2.1 Programming sequence in free mode

Step 1

The first programme entered by the Thermo 500 after pressing the free mode key  is the “Free Design” programme. It is a standard programme which can be changed and executed, but which is never saved after turning off the device.

Use the turning knob to select a free programme position or a previously saved programme. The free programmes are numbered from 01 to 49.

Caution!

In Free mode, indication and standard mode the Thermo 500 automatically detects the type of Thermoplane being connected. You can choose which type of Thermoplane you wish and to which channel you want to connect the Thermoplane.

If you opt for the free mode, which has a higher average power than 32 Watt and you connect a Thermoplane 80, the Thermo 500 will automatically reduce the average power to 32 Watt. This will reduce the peak power. In that case a flashing Thermoplane 140/80 symbol indicates that you have connected the Thermoplane 80.

The average power of the Thermoplane 80 cannot be set higher than 32 Watt. The Thermo 500 will always check this. If, for example, you set the frequency higher, so that the average power of 32 Watt is exceeded, the peak power will automatically be reduced.

The formula that is used to calculate this is:

$$\text{Frequency} \times \text{pulse time} \times \text{peak power} \leq 32 \text{ Watt}$$

Example:

4	0	0	μ	s				4	0	0	H	z				6	5	μ	s
2	0	0	W					2	0	m	i	n				1	0	W	

According to the formula the average power is precisely 32 Watt.

The pulse frequency is then increased to 500 Hz and the peak power is automatically adjusted to 160 Watt so that the average power remains 32 Watt.

4	0	0	μ	s				5	0	0	H	z				6	5	μ	s
1	6	0	W					2	0	m	i	n				1	0	W	

Caution

A parameter is saved automatically when changed.

Pay attention not to overwrite an existing programme unintentionally.

The Thermo 500 automatically detects the presence of an Thermoplane and its characteristics.

Step 2: Programme name

The name of the standard programme is “FREE N”, N stands for the number of the following free memory place . You can overwrite or keep this name. If you want to keep the name, push the SELECT > key  till you reach the following parameter, without changing the letters by means of the turning knob.

To overwrite the name, select the letters to be changed, by pressing SELECT  and change them by means of the turning knob (blank, A...Z, 0...9, _ , - , and so on.):

						F	R	E	E		0	1					
0	.	0	0	W						2	0	m	i	n			

		P	_			F	R	E	E		0	1					
0	.	0	0	W						2	0	m	i	n			

		P	R	_		F	R	E	E		0	1					
0	.	0	0	W						2	0	m	i	n			

And so on till the desired name, for example “PROGRAM 4” is ready.

Step 3

Push the SELECT > key  to select the parameter of the treatment time. The pre-selection is 20 minutes.

		P	R	O	G	R	A	M		4							
0	.	0	0	W						2	0	m	i	n			

The treatment time parameter can be set in steps of 1 minute till a maximum value of 60 minutes.

In the START window no other parameters can be set. In this way the programme place is again selected by pressing SELECT  again.

Step 4

Push the free mode key  to display the parameter window. When a Thermopode of the right channel has also been connected, these parameters are also displayed on the screen.

First the pulse time of the left Thermopode is selected.

Attention! All the Free mode programs have already the minimum parameter setting.

	6	5	μ	s				2	5	H	z				6	5	μ	s
	1	0	W					2	0	m	i	n			1	0	W	

The pulse time of the left channel can be set from 65 μs to 400 μs in steps of 5 μs.

Step 5

Push the SELECT > key  to select the parameter pulse frequency.

	6	5	μ	s				2	5	H	z				6	5	μ	s
	1	0	W					2	0	m	i	n			1	0	W	

The pulse frequency can be set

- from 25 Hz to 875 Hz in one channel therapy **and** in connecting the Thermoplude 140 or
- from 25 Hz to 500 Hz in two channel therapy **and** in connecting two Thermoplodes 140
- from 25 Hz to 875 Hz in connecting a Thermoplude 80
- from 25 Hz to 500 Hz in connecting two Thermoplodes 80
- from 25 Hz to 500 Hz in connecting a Thermoplude 140 and a Thermoplude 80.

Setting in steps of 5 Hz. **This parameter pulse frequency applies to both channels!**

Step 6

Push SELECT  2x to select the peak power parameter of the left channel.

	6	5	μ	s				2	5	H	z				6	5	μ	s
	1	0	W					2	0	m	i	n			1	0	W	

By means of the turning knob you can set the peak power from 10 W to 200 W in steps of 5 W. As long as the start/stop LED is off, you cannot start the programme. In any case the peak power has to be set in order to have the start/stop LED lightened up.

Caution!

If a Thermoplude 80 is connected, the peak power to be set will be limited if the average power of 32 Watt is exceeded. (This occurs for high parameters for the pulse frequency and pulse time.)

Step 7

Push the SELECT > key  possibly to select the treatment time parameter, which in this window can also be set in a previously started programme. The time parameter can be set in steps of 1 minute, till 60 minutes maximally.

All the parameters set in stand-by mode are automatically saved.

Programming illogical parameters is impossible. As long as the start/stop LED is off, you cannot start the programme.

Step 8

The programme can be started by the start/stop key . The START window is shown on the display. After the start the average power first equals 0.00 W.

		P	R	O	G	R	A	M	4							
0	.	0	0	W				2	0	:	0	0				

In practice the treatment only starts by turning on the power knob. Only at that moment the treatment time starts to run. The turning knob is automatically active for the left channel.

When the running (active) programme disposes of a right channel function and in case the right channel Thermoplude is connected, the average power of the right channel can be selected by pressing the SELECT > key . Following the selection of the average power parameter of the right channel, you can set it by means of the turning knob.

Step 9

When you want to change a parameter during the treatment, you have to push the

free mode key  again to display the parameter window. To illustrate this, the parameter pulse time of the left channel is selected:

4	0	0	μ	s			4	0	0	H	z				6	5	μ	s	
	2	0	0	W				1	9	:	3	0				1	0	W	

The parameter can be set by means of the turning knob. The setting directly affects the average power and the display of Prms in the START window.

Push the SELECT > key  when you want to set the next parameter.

During the treatment (meaning with an increased power) the changed/modified parameters are not saved !

At any time the START window can be called by pressing the free mode key .

At the end of the treatment time and by pressing the start/stop key  the treatment is terminated and the parameters are set back to the initial value of the programme.
At the end of the treatment time you will also hear a buzzing sound of about 10 seconds.

4.3.3 Execution of the programme

The desired programme is called by means of the turning knob.

		P	R	O	G	R	A	M		1	6								
3	2	.	0	W				1	0	:	0	0		1	6	.	0	W	

Push the start/stop key  to start the programme. You can continue working as described in step 10 of paragraph 4.3.2.1.

5 CONTROL OF THE ELECTRODE CONNECTION

The Thermo 500 disposes of an on-line control of the connection and type of Thermopode.

In Free mode, indication and standard mode the Thermo 500 automatically detects the type of Thermopode being connected. You can choose which type of Thermopode you wish and to which channel you want to connect the Thermopode.

If you opt for the free mode or indication which has a higher average power than 32 Watt and you connect a Thermopode 80, the Thermo 500 will automatically reduce the average power to 32 Watt. This will reduce the peak power and/or pulse frequency.

A flashing Thermopode 140/80 symbol indicates that you have connected the Thermopode 80.

The average power of the Thermopode 80 cannot be set higher than 32 Watt. The Thermo 500 will always check this. If, for example, you set the frequency higher (max. 500 Hz), the peak power will automatically be reduced.

The formula that is used to calculate this is:

Frequency x pulse time x peak power $\leq$ 32 Watt

Example :

4	0	0	μ	s			4	0	0	H	z					6	5	μ	s
2	0	0	W					2	0	m	i	n				1	0	W	

If the pulse frequency is then increased to 500 Hz then the peak power will automatically be adjusted to 160 Watt

4	0	0	μ	s			5	0	0	H	z					6	5	μ	s
1	6	0	W					2	0	m	i	n				1	0	W	

6 TREATMENT

6.1 Preparing the patients

To achieve an optimal way of treatment, the patients and the body part to be treated should best be relaxed completely. Lying or sitting comfortably is essential. Patients cannot be treated in metal chairs or couches. For safety reasons it is recommended that hearing-aids, watches, rings, necklaces, bracelets and other metal objects be removed before the treatment. The same principle applies to people who are in the vicinity of the appliance or near the Thermoplodes and cables, for example the medical staff.

The body parts concerned are treated naked. Clothes made of synthetic fabrics have to be removed because of the low absorbing capacity. Moist skin, for example in skin creases, can cause local overheating. To avoid this and from a hygienic point of view, it is recommended to put a dry layer of cellulose or a thin terry cloth between the Thermoplude and the treated body part. The patient cannot touch metal objects during the treatment, neither the appliance itself, nor the metal parts of connections.

6.2 Thermal effect of the coil-field method

The coil-field electrodes cause Foucault currents as a result of their high-frequency magnetic field, converted into heat in the tissue. The currents are even higher, when the electric conduction of the tissue is higher (tissue with a good blood circulation, for example muscular tissue and internal organs). To reach the deeper located tissue, the coil-field electrodes of the Thermo 500 have an electrostatic fence reducing the electric field of the coil-field electrodes by which the coil-field electrodes which impede the thermal effect on the epidermis tissue.

6.3 Indications

Specific indications for pulsing short-wave therapy are:

Post-traumatic lesions, such as:

- distortion
- contusion
- rupture
- fracture
- haematoma

Post-operative lesions

Inflammation, such as :

- arthritis (in a non-active stage)
- bursitis, possibly with calcification
- sinusitis
- tendinitis
- wound healing

Degenerative lesions, such as :

- arthrosis
- rheumatoid arthritis

Lesions of the peripheral nerves, such as :

- pareses
- neuritoids, neuralgias
- herpes zoster
- neuropathy

6.3.1 Contra-indications

General contra-indications:

- high fever
- serious cardiac disorders
- psychological problems (aversion of the patient, fear)
- general : cancer, malignant tumours (metastasis)
- general : tuberculosis.
- in principle : reduced heat sensation in the affected area

Implanted metals

It is not allowed to treat patients with pacemakers with high-frequency heat therapy. Body parts with metal parts (pins, shell shrapnel or similar things) cannot be treated.

Pregnant women, threatened organs

The treatment of pregnant women in the abdomen is not advised.

Other contra-indications are :

- . growing zones of the bone/os
- . malignant tumours
- . tuberculosis
- . arterial haemorrhage disorders of stage III and IV
- . varicose veins
- . general tendency to bleed

In organs with a minor vasculature and blood circulation (eyes, testicles) administer the dose carefully!
Also the thermal dose with patients with a disordered heat sensation must be avoided.

7 MALFUNCTIONS, GUARANTEE, RELIABILITY, CUSTOMER CARE

7.1 Malfunctions

The Thermo 500 independently tests itself. If functions do not work anymore, so that the parameters are situated outside the allowed deviations, the mention “Call Customer Service” appears on the display. This mention can also appear on the base of other (external) factors. After switching the appliance on and off, it is ready for use again.

Faults, which are caused by a wrong use of the appliance, are clearly indicated on the display.

7.2 Guarantee and reliability

7.2.1 Terms of guarantee

The guarantee does not apply to repairs of faults or defects originated by a negligent use of the appliance, by a wrong interpretation of the faults or a non-correspondence to the instructions of this manual, by not observing the instruction manual, or by accidents, caused by maintenance or repair by people not authorised by the company GymnaUniphy.

7.2.2 Liability of the manufacturer

After a period of 10 years following the first operation of the appliance (or its accessories) the manufacturer cannot be held liable for defects to the appliance, to its accessories, or for possible ensuing damage.

Neither can the manufacturer be held liable for ensuing injury or harm to therapists, patients or the used accommodation after for example a wrong diagnosis, non-professional utilisation of the appliance or its accessories, wrong interpretation or non-observance of the instruction manual, after a bad maintenance of the appliance or in case a maintenance or repair have been carried out by persons not authorised by the manufacturer.

The manufacturer, nor the supplier, can be held responsible by any means for the transmission of infections through catheters/probes or electrodes.

7.3 Customer care

Your supplier is only responsible for effective operation of the appliance when:

- Repairs are carried out by authorised persons;
- the electrical installations in the area involved comply with the applicable legal guidelines;
- all repairs, modifications, extensions or settings are carried out by authorised persons according to the guidelines of this instruction manual;
- the appliance is used for the purpose for which it has been designed;
- maintenance has been carried out regularly and in the described manner;
- the legal guidelines for using the appliance have been observed.

Apart from the fuses, no components in the appliance can be replaced by the user. In the event of improper use or lack of maintenance according to the guidelines, GymnaUniphy and its representatives are released of all liability of any ensuing damage, injury, breakdowns and malfunctions.

Your supplier takes care of all the interventions as to customer service and guarantee.

Your supplier's terms and conditions of delivery are applicable.

The guarantee lapses if the appliance is not used in accordance with the instructions above. The expected lifespan of the appliance is 10 years.

8 MAINTENANCE AND CLEANING

For maintenance and cleaning to the Thermo 500 first switch off the appliance and remove the plug from the wall socket.

8.1 Cleaning and disinfection

For the normal, daily maintenance of the appliance, the appliance and its accessories need not to be opened. If necessary the outside can be cleaned by means of a dry or slightly moistened cloth.

We strictly advise against the use of solvents, considering they can damage the used materials. Persistent stains can be removed with a non-caustic soap substance.

Anyway prevent any liquid from entering the appliance or its accessories! All parts must be dried carefully.

To clean and disinfect the appliance and its accessories (except for felt parts) we advise disinfectants on sale everywhere, according to the instructions. For erasing and sprinkling disinfection, only disinfectants with aldehyde, alcohol or quaternary ammonia compounds can be used, for protective reasons, observing the guidelines of the manufacturer regarding dilution and effect.

For users in Germany we only advise the use of disinfectants appearing on the latest DGHM list (DGHM: Deutsche Gesellschaft für Hygiene und Mikrobiologie). This list also mentions the substance basis of the disinfectants. Because of possible damage to the materials, preparations on the base of halogen secreting compounds cannot be used, nor strong organic acids, compounds giving off oxygen, solvents, petrol and so on...

To guarantee the quality and safety, unauthorised persons cannot open the appliance or its accessories by any means. Opening the appliance or its accessories (for maintenance or repair) is preserved to qualified personnel authorised by the company GymnaUniphy.

9 TECHNICAL SPECIFICATIONS

Current connection	230 V $\pm$ 10 %: 50/60 Hz 115 V $\pm$ 10 %; 50/60 Hz
Electric fuse/external	10 A at 230 V 16 A at 115 V
Fuses, externally accessible	6,3 AF at 230 V 6,3 AT at 115 V
Insulation class	I, Type BF (with connected electrodes)
Safety degree	IP 20
Power absorption max.	Approximately 450 VA
Power in stand-by	Approximately 100 VA
Frequency	27,12 MHz $\pm$ 0,6 %
Peak power P _{pulse}	200 W at 50 Ohm
Average power/power P_{rms}	
- Per channel, Thermoplude 140 in two-channel operation	Max. 40 W
- 1 channel, Thermoplude 140	Max. 70 W
- Per channel, Thermoplude 80 in two-channel operation	Max. 32 W
- Thermoplude 80 and Thermoplude 140 in two-channel operation	Thermoplude 80 : Max. 32 W Thermoplude 140 : Max. 40 W
Pulse time	65 – 400 μ s
Pulse frequency	25 - 500 Hz in two channel operation 25 - 875 Hz in one channel operation
Treatment time	0 – 60 minutes
Dimensions (w x d x h)	371 mm x 356 mm x 930 mm
Weight	Approximately 43 kg (excl. Thermoplude)
Exits	2 coaxial exits (protected) for Thermoplodes
Standards	EN 60 601-1: 1990+A1+A2+A13 EN 60 601-1-2: 1993 EN 60 601-2-3: 1993+A1
MPG-class	IIa
EC sign	complying with the directive MDD: 93/42/EEC

10 SURVEY OF THE ACCESSORIES

10.1 Standard accessories

<u>Number</u>	<u>Description</u>
1	instruction manual Thermo 500
1	coil-field electrode “Thermoplode” 140, 14 cm diameter, excl. cable
1	electrode cable for the “Thermoplode”
1	electrode arm with a 10 mm open-end wrench to set the arm
1	Swiss adaptor

10.2 Optional accessories

Description

- Coil-field electrode “Thermoplode 80”, diameter 8 cm, excl. cable
- Electrode cable for the “Thermoplode”
- Coil-field electrode “Thermoplode 140” diameter 14 cm, excl. cable
- Electrode arm with a 10 mm open-end wrench to set the arm

We are always concerned to keep our appliances technically up-to-date. Therefore we entitle ourselves to carry out modifications or adaptations whenever necessary.

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